

SAFETY DATA SHEET

ANALYTICAL REFERENCE STANDARD

Date created: November 28, 2017

Date reviewed: N/A

SECTION 1 - Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name: alpha-Hydroxyalprazolam-D5 (0.1 mg base/1 ml in methanol)
Article number: ALP-923-FB-0.1LM

1.2 Relevant identified uses of the substance or mixture and uses advised against

Reference material for analytical purposes only.

1.3 Details of the supplier of the safety data sheet

Lipomed AG
Fabrikmattenweg 4
CH-4144 Arlesheim

Tel: +41 (0) 61 702 02 00
Fax: +41 (0) 61 702 02 20
eMail: lipomed@lipomed.com

1.4 Emergency phone number

Tel: +41 (0) 44 321 34 33

SECTION 2 - Hazard identification

2.1 Classification of the substance or mixture

Classification according to European Regulation (EC) 1272/2008 (EU-GHS/CLP)

Flammable liquids (Cat. 2)	H225: Highly flammable liquid and vapour.
Acute toxicity oral (Cat. 3)	H301: Toxic if swallowed.
Acute toxicity dermal (Cat. 3)	H311: Toxic in contact with skin.
Acute toxicity inhalation (Cat. 3)	H331: Toxic if inhaled.
Specific target organ toxicity (Cat. 1) Single exposure	H370: Causes damage to organs.

2.2 Label elements

Labelling according to European Regulation (EC) 1272/2008 (CLP)

Hazard Pictograms:



GHS02



GHS06



GHS08

Signal word: Danger.

Hazard Statements:

H225 Highly flammable liquid and vapour.
H301 Toxic if swallowed.

- H311 Toxic in contact with skin.
- H331 Toxic if inhaled.
- H370 Causes damage to organs.

Precautionary Statements:

Prevention

- P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking.
- P260 Do not breathe dust/fume/gas/mist/vapors/spray.
- P264 Wash skin thoroughly after handling.
- P270 Do not eat, drink or smoke when using this product.
- P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response

- P301+P310 IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.
- P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
- P303+P361+P353 IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.
- P307+P311 If Exposed: Call a POISON CENTER or doctor/physician.
- P361 Remove/ Take off immediately all contaminated clothing.
- P405 Store locked up.

2.3 Others hazards

None.

SECTION 3 - Composition/information on ingredients

3.1 Substances

Not applicable.

3.2 Mixtures

alpha-Hydroxyalprazolam-D5 (0.1 mg base/1 ml in methanol)

Product Name	Methanol	alpha-Hydroxyalprazolam-D5
Formula	CH4O	C17H8D5CIN4O
Molecular Weight	32.04	329.79
CAS Number	67-56-1	136765-24-7
EC Number	200-659-6	802-922-6
Percent	99.99	0.01

SECTION 4 - First aid measures

4.1 Description of first aid measures

FOLLOWING ANY ROUTE OF EXPOSURE GET MEDICAL ATTENTION IMMEDIATELY.



Lipomed AG is **ISO 9001** certified and **ISO 17025/34** accredited.

SERIOUS TOXICITY IS PRECEDED BY VOMITING IN MOST CASES OF ORAL INGESTION.

Inhalation:

If inhaled, remove to fresh air. If breathing is labored or with coughing, give 100% supplemental oxygen. If not breathing, begin artificial respiration. DO NOT GIVE MOUTH-TO-MOUTH RESUSCITATION.

Ingestion:

If swallowed, get medical attention immediately; do not induce vomiting. Never give anything by mouth to an unconscious person. If not breathing, begin artificial respiration. Call a physician immediately.

Skin Contact:

Immediately flush skin with plenty of water for at least 15 minutes while removing contaminated clothing and shoes. Call a physician immediately.

Eye Contact:

Immediately flush eyes with plenty of water for at least 15 minutes, lifting lower and upper eyelids occasionally. Call a physician immediately.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in section 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available.

SECTION 5 - Firefighting measures

5.1 Extinguishing Media

Suitable: Water spray; carbon dioxide; dry chemical powder; appropriate foam.

5.2 Special hazards arising from the substance or mixture

Specific Hazard: Flammable liquid, emit toxic fumes under fire conditions.

Explosion Hazard: Vapor may travel considerable distance to source of ignition and flash back.

5.3 Advice for firefighters

Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

SECTION 6 - Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear protective equipment. Avoid breathing vapors. Shut off all sources of ignition. Keep unprotected persons away.

6.2 Environmental precautions

Dilute with plenty of water. Do not allow material to enter drains and water sources.



6.3 Methods and material for containment and cleaning up

To clean up, cover with dry-lime, sand or soda ash. Place in covered containers using non-sparking tools and transport outdoors. Wash spill site with appropriate cleaning agent to removal residual traces of spilled material.

6.4 Reference to other sections

See section 13 for disposal.

SECTION 7 - Handling and storage

7.1 Precautions for safe handling

Do not breathe vapor. Do not get in eyes, on skin, on clothing. Avoid prolonged and repeated exposure. Ensure good ventilation at workplace. Do not smoke.

7.2 Conditions for safe storage, including any incompatibilities

Keep tightly closed. Store in a cold and dry place. Protect from light. Keep away from heat, sparks and open flame. Refer to the manufacturer's certificate for storage conditions.

7.3 Specific end use(s)

No specific information.

SECTION 8 - Exposure controls/personal protection

8.1 Control parameters

Methanol: requires monitoring at working place.

Exposure limits: Methanol

REL: Short term value: 250 ppm
Long term value: 200 ppm

TLV: Short term value: 250 ppm
Long term value: 200 ppm

8.2 Exposure controls

Skin Protection:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.

The protective gloves to be used must comply with the specifications of EC Directive 89/686/EEC and the related standard EN374.

Wear impervious protective clothing, including boots, gloves, lab coat, apron or coveralls, as appropriate, to prevent skin contact. Dispose contaminated material according good laboratory practices.

Eye/face Protection:

Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Personal Respirators (NIOSH Approved):

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multipurpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls.

If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

SECTION 9 - Physical and chemical properties

9.1 Information on basic physical and chemical properties

Methanol:

Property	Value
Physical state:	
Colour	Colorless
Form	Liquid
Odour	No data available
Melting point/freezing point	- 98 °C/ No data available
Boiling point or initial boiling point and boiling range	64.6 °C (64 – 65°C) at 1.013 hPa (0.760 mmHg)
Flammability	No data available
Lower and upper explosion limit	Lower explosion limit: 6 % (V) Upper explosion limit: 36 % (V)
Flash point	11 °C
Auto ignition temperature	No data available
Decomposition temperature	No data available
pH	No data available
Kinematic viscosity	No data available
Partition coefficient n-octanol/water	No data available
Vapor pressure	13.02 kPa at 20°C
Density and/or relative density	0.791 g/cm ³ at 20°C
Relative vapor density	No data available

9.2 Other information

No data available.

SECTION 10 - Stability and reactivity

10.1 Reactivity

No data available.

10.2 Chemical stability

Stable under recommended storage conditions.



10.3 Possibility of hazardous reactions

No data available.

10.4 Conditions to avoid

Sources of ignition: Heat, flames.

10.5 Incompatible materials

Strong acids, strong oxidizing agents, alkali metals, strong reducing agents.

10.6 Hazardous decomposition products:

Possible formation of carbon oxides during heat and in case of fire.

SECTION 11 - Toxicological information

11.1 Information on toxicological effects

Acute toxicity (methanol):

oral; rat: LD50: 5628 mg/kg

inhalation; rat: LC50: 64000 ppm/4h

intraperitoneal; rat: LD50: 7529 mg/kg

oral; mouse: LD50: 7300 mg/kg

intraperitoneal; mouse: LD50: 10765 mg/kg

Skin corrosion/irritation

May cause skin irritation.

Serious eye damage/irritation

May cause eye damage/irritation.

Respiration or skin sensitization

May cause respiration or skin sensitization.

Germ cell mutagenicity

No data available.

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

Reproductive toxicity

No data available.

Specific target organ toxicity (STOT) – single exposure

No data available.

Specific target organ toxicity (STOT) – repeated exposure

No data available.

Aspiration hazard

No data available.

Additional information

RTECS (methanol): PC1400000

Formulation not fully tested yet.

Methanol may be fatal or cause blindness is swallowed.

Effects due to Ingestion may include: Headache, vomiting, Nausea, dizziness.

SECTION 12 - Ecological information

12.1 Toxicity

Methanol

Fish: LC50: 24000 mg/l/96h.

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

Not applicable.

12.6 Other adverse effects

No data available.

SECTION 13 - Disposal considerations

13.1 Waste treatment methods

Contact a licensed professional waste disposal service to dispose of this material. Burn in a chemical incinerator equipped with an afterburner and scrubber. Observe all federal, state, and local environmental regulations but exert extra care in igniting as this material is highly flammable. Observe all federal, state, and local environmental regulations.

SECTION 14 – Transport information

14.1 UN number

RID/ADR: 1230

IMDG: 1230

IATA: 1230

14.2 UN proper shipping name

RID/ADR: Methanol (Mixture)

IMDG: Methanol (Mixture)

IATA: Methanol (Mixture)

14.3 Transport hazard class(es)

RID/ADR: 3

IMDG: 3

IATA: 3

Sub. Risk: 6.1

Sub. Risk: 6.1

Sub. Risk: 6.1

14.4 Packaging group

RID/ADR: II

IMDG: II

IATA: II

14.5 Environmental hazards

RID/ADR: No

IMDG: No marine pollutant

IATA: No



Lipomed AG is **ISO 9001** certified and **ISO 17025/34** accredited.

14.6 Special precautions for user

RID/ADR: Small quantities
exempted

IMDG: Small quantities
exempted

IATA: Small quantities
exempted

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Not applicable.

SECTION 15 - Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This SDS complies with Regulation: No.: 1907/2006.

15.2 Chemical safety assessment

Chemical Safety Assessment has not been carried out.

SECTION 16 - Other information

This information is given without any warranty or representation. It is believed to be correct but does not claim to be all inclusive and shall be used only as a guide. Lipomed AG shall not be held liable for any damage resulting from handling or contact with the above products. It is offered solely for your consideration, investigation and verification.

DISCLAIMER

For R & D only. Not for drug, household or other uses.



Lipomed AG is **ISO 9001** certified and **ISO 17025/34** accredited.