




» STEROIDS

The use and abuse of performance-enhance substances, e.g. anabolic steroids, continues to plague competitive sports. The detection of synthetic steroids is a continuing challenge for analytical laboratories especially because of the availability of numerous of such compounds.

As for other drugs and/or drugs of abuse, there is a major interest to measure steroids in biological samples for sport but also for clinical and forensic cases.

To support laboratories involved in the steroid testing, LIPOMED now supply high quality reference standards for the steroids of interest available as ready-to-use calibrated solution and/or as standard pack sizes of 10 mg, 50 mg and 100 mg.

For more information, please check our website www.lipomed.com.




Androstenedione

Synonym : 4-Androstene-3,17-dione

$C_{19}H_{26}O_2$ M.W.: 286.41

AND-1489-1LM	1.0 mg/1 ml	calibrated solution in methanol
AND-1489-10	10 mg	powder
AND-1489-50	50 mg	powder
AND-1489-100	100 mg	powder




Boldenone

Synonym : 1-Dehydrotestosterone

$C_{19}H_{26}O_2$ M.W.: 286.42


BOL-1581-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
BOL-1581-10	10 mg	powder
BOL-1581-50	50 mg	powder
BOL-1581-100	100 mg	powder



Chlorodehydromethyltestosterone

$C_{20}H_{27}ClO_2$ M.W.: 334.88

CDM-1813-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
CDM-1813-10	10 mg	powder
CDM-1813-50	50 mg	powder
CDM-1813-100	100 mg	powder




Clostebol acetate

Synonym : 4-chlorotestosterone 17 β -acetate (4-CLTA)

$C_{21}H_{29}ClO_3$ M.W.: 364.91

CLS-1816-10	10 mg	powder
CLS-1816-50	50 mg	powder
CLS-1816-100	100 mg	powder




Corticosterone (shipment under ice for solution)

Synonym : 17-Deoxycortisol

$C_{21}H_{30}O_4$ M.W.: 346.47

COR-1628-1LM	1.0 mg/1 ml	calibrated solution in methanol
COR-1628-10	10 mg	powder
COR-1628-50	50 mg	powder
COR-1628-100	100 mg	powder




Cortisol (shipment under ice for solution)

Synonym : Hydrocortisone

$C_{21}H_{30}O_5$ M.W.: 362.47


COR-1630-1LM	1.0 mg/1 ml	calibrated solution in methanol
COR-1630-10	10 mg	powder
COR-1630-50	50 mg	powder
COR-1630-100	100 mg	powder



Cortisone (shipment under ice for solution)

$C_{21}H_{28}O_5$ M.W.: 360.45

COR-1633-1LM	1.0 mg/1 ml	calibrated solution in methanol
COR-1633-10	10 mg	powder
COR-1633-50	50 mg	powder
COR-1633-100	100 mg	powder




Dehydroepiandrosterone

Synonym : DHEA

$C_{19}H_{28}O_2$ M.W.: 288.42

DHE-1603-1LM	1.0 mg/1 ml	calibrated solution in methanol
DHE-1603-10	10 mg	powder
DHE-1603-50	50 mg	powder
DHE-1603-100	100 mg	powder




11-Deoxycortisol

Synonym : Cortexolone; Cortodoxone

$C_{21}H_{30}O_4$ M.W.: 346.47

COR-1627-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
COR-1627-10	10 mg	powder
COR-1627-50	50 mg	powder
COR-1627-100	100 mg	powder




5 alpha-Dihydrotestosterone (shipment under ice for solution)

Synonym : DHT, Androstanolone

$C_{19}H_{30}O_2$ M.W.: 290.45

DHT-1685-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
DHT-1685-10	10 mg	powder
DHT-1685-50	50 mg	powder
DHT-1685-100	100 mg	powder




Drostanolone

Synonym : Masteron

$C_{20}H_{32}O_2$ M.W.: 304.47


DRO-1710-10	10 mg	powder
DRO-1710-50	50 mg	powder
DRO-1710-100	100 mg	powder



17 beta-Estradiol


$C_{18}H_{24}O_2$ M.W.: 272.39

EST-1684-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
EST-1684-10	10 mg	powder
EST-1684-50	50 mg	powder
EST-1684-100	100 mg	powder




Fluoxymesterone

$C_{20}H_{29}FO_3$	M.W.: 336.44	
FMT-1814-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
FMT-1814-10	10 mg	powder
FMT-1814-50	50 mg	powder
FMT-1814-100	100 mg	powder




17α-Hydroxyprogesterone

$C_{21}H_{30}O_3$	M.W.: 330.47	
PRO-1629-1LM	1.0 mg/1 ml	calibrated solution in methanol
PRO-1629-10	10 mg	powder
PRO-1629-50	50 mg	powder
PRO-1629-100	100 mg	powder




6-β-Hydroxytestosterone

$C_{19}H_{28}O_3$	M.W.: 304.43	
TES-843-1LM	1.0 mg/1 ml	calibrated solution in methanol
TES-843-10	10 mg	powder
TES-843-50	50 mg	powder
TES-843-100	100 mg	powder



Metenolone acetate

Synonym : Methenolone acetate		
$C_{22}H_{32}O_3$	M.W.: 344.49	
MTN-1817-1LI	1.0 mg/1 ml	calibrated solution in isopropanol
MTN-1817-10	10 mg	powder
MTN-1817-50	50 mg	powder
MTN-1817-100	100 mg	powder




Metenolone enanthate

Synonym : Methenolone enanthate

$C_{27}H_{42}O_3$ M.W.: 414.62


MTN-1818-1LI	1.0 mg/1 ml	calibrated solution in isopropanol
MTN-1818-10	10 mg	powder
MTN-1818-50	50 mg	powder
MTN-1818-100	100 mg	powder



Methandienone

$C_{20}H_{28}O_2$ M.W.: 300.44


TES-1249-1LM	1.0 mg/1 ml	calibrated solution in methanol
TES-1249-10	10 mg	powder
TES-1249-50	50 mg	powder
TES-1249-100	100 mg	powder



17alpha-Methyltestosterone

$C_{20}H_{30}O_2$ M.W.: 302.46

TES-1810-1LM	1.0 mg/1 ml	calibrated solution in methanol
TES-1810-10	10 mg	powder
TES-1810-50	50 mg	powder
TES-1810-100	100 mg	powder




Nandrolone

Synonym : 19-Nortestosterone

$C_{18}H_{26}O_2$ M.W.: 274.40

TES-889-1LM	1.0 mg/1 ml	calibrated solution in methanol
TES-889-10	10 mg	powder
TES-889-50	50 mg	powder
TES-889-100	100 mg	powder




Oxymetholone (shipment under ice for solution)

Synonym : Anasteron

$C_{21}H_{32}O_3$ M.W.: 332.48


OXY-1678-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
OXY-1678-10	10 mg	powder
OXY-1678-50	50 mg	powder
OXY-1678-100	100 mg	powder



Progesterone

$C_{21}H_{30}O_2$ M.W.: 314.47


PRO-1683-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
PRO-1683-10	10 mg	powder
PRO-1683-50	50 mg	powder
PRO-1683-100	100 mg	powder



Stanozolol

$C_{21}H_{32}N_2O$ M.W.: 328.50

STA-1522-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
STA-1522-10	10 mg	powder
STA-1522-50	50 mg	powder
STA-1522-100	100 mg	powder




Testosterone

Synonym : trans-Testosterone


$C_{19}H_{28}O_2$ M.W.: 288.42

TES-842-1LM	1.0 mg/1 ml	calibrated solution in methanol
TES-842-10	10 mg	powder
TES-842-50	50 mg	powder
TES-842-100	100 mg	powder




Testosterone-D₃

$C_{19}H_{25}D_3O_2$	M.W.: 291.44	
TES-1149-0.1LM	0.1 mg/1 ml	calibrated solution in methanol
TES-1149-1LM	1.0 mg/1 ml	calibrated solution in methanol
TES-1149-10	10 mg	powder
TES-1149-50	50 mg	powder
TES-1149-100	100 mg	powder




Testosterone 17-heptanoate

Synonym : Testosterone 17-enanthatate		
$C_{26}H_{40}O_3$	M.W.: 400.59	
TES-1242-10	10 mg	powder
TES-1242-50	50 mg	powder
TES-1242-100	100 mg	powder




Testosterone 17-isocaproate

$C_{25}H_{38}O_3$	M.W.: 386.57	
TES-1243-10	10 mg	powder
TES-1243-50	50 mg	powder
TES-1243-100	100 mg	powder



Testosterone 17-phenylpropionate


$C_{28}H_{36}O_3$	M.W.: 420.56	
TES-1241-10	10 mg	powder
TES-1241-50	50 mg	powder
TES-1241-100	100 mg	powder



Testosterone 17-propionate

$C_{22}H_{32}O_3$ M.W.: 344.49


TES-1240-10	10 mg	powder
TES-1240-50	50 mg	powder
TES-1240-100	100 mg	powder



Testosterone 17-undecanoate

$C_{30}H_{48}O_3$ M.W.: 456.70

TES-1253-10	10 mg	powder
TES-1253-50	50 mg	powder
TES-1253-100	100 mg	powder




Trenbolone

Synonym : 17 beta-Trenbolone

$C_{18}H_{22}O_2$ M.W.: 270.37

TBL-1529-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
TBL-1529-10	10 mg	powder
TBL-1529-50	50 mg	powder
TBL-1529-100	100 mg	powder



Trenbolone acetate

Synonym : 17 beta-Trenbolone acetate

$C_{20}H_{24}O_3$ M.W.: 312.41

TBL-1528-1LI	1.0 mg/1 ml	calibrated solution in isopropanol
TBL-1528-10	10 mg	powder
TBL-1528-50	50 mg	powder
TBL-1528-100	100 mg	powder

Ordering & other information

LIPOMED is a private Swiss Health Care Company. We are one of the leading companies for Analytical Reference Standards worldwide and specialized in the research, development and manufacturing of ethical pharmaceuticals and products for the treatment of rare diseases.

Since 1993 LIPOMED is dedicated to providing products and services of outstanding quality and expertise. We have our own GMP/GDP and ISO 9001 certified pharmaceutical production facilities in Arlesheim, near Basel. We are also ISO/IEC 17025 accredited for testing of Analytical Reference Standards and ISO 17034 accredited for the production of Reference Material.

Our innovative products are marketed and distributed with offices and partners in more than 50 countries around the world.

LIPOMED's Analytical Reference Standards currently include the following classes of substances, their deuterated analogues and metabolites, such as Glucuronides and Sulfates:

- » Alcohol & Ethanol Biomarkers
- » Amphetamines & Sympathomimetic amines
- » Antidepressants
- » Antipsychotics
- » Barbiturates
- » Benzodiazepines
- » Cannabinoids
- » Cathinones
- » Cocaines
- » Custom Mixtures
- » Hallucinogens
- » Opioids
- » Other CNS Depressants
- » Steroids
- » Vitamins

For more information about LIPOMED's Analytical Reference Standards, please visit our website www.lipomed.com.

SWITZERLAND (HEADQUARTERS)

LIPOMED AG
Fabrikmattenweg 4
4144 Arlesheim
Switzerland
Tel +41 61 702 02 00
Fax +41 61 702 02 20
lipomed@lipomed.com

GERMANY

LIPOMED GmbH
Hegenheimer Str. 2
79576 Weil am Rhein
Germany
Tel +49 7621 1693 473
Fax +49 7621 1693 474
lipomed@lipomed.com

USA, CANADA, BERMUDA

LIPOMED Inc
150 Cambridgepark Drive
Cambridge, MA 02140
USA
Tel +1 617 577 7222
Fax +1 617 577 1776
ussales@lipomed.com



LIPOMED AG, Switzerland is **GMP/GDP** and **ISO 9001** certified for production and distribution of Pharmaceuticals and **ISO/IEC 17025** (ANAB certificate no. AT-1760) accredited for testing of Analytical Reference Standards and **ISO 17034** (ANAB certificate no. AR-1761) accredited for Reference Material production. All products distributed by LIPOMED affiliates are manufactured by LIPOMED AG and distributed under protocols of LIPOMED AG.

